







Test report



At-home test



IBS Test

Lab test

Stool

Name: **Sample Report** Date of test: **09/28/2023** Analysis-ID: **DUMMY-41**

Table of contents



















Determination of aerobic bacteria	page 7
Determination of anaerobic bacteria	page 8
Mycological stool examination	page 10
pH-value	page 11
Digestive residues	page 11
Detection of indigestion	page 12
Detection of malabsorption	page 13
Mucosal immunity	page 14

Your test results









How to interpret your results

If the score on the scale is within the green area, your value is good. If the score is within the yellow, orange or red area, it indicates a deviation. Some parameters can not be set too high or too low, in which case the scale starts or ends on the green area, this is completely correct. For an explanation of the various parameters please see Part 2 entitled Extended information.







Determination of aerobic bacteria

Name	Your value	Unit	Reference value	Scale
Escherichia coli	 8,0 x 10 ⁵	CFU/g faeces	10 ⁶ - 10 ⁷	
Escherichia coli Biovare	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Proteus spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Klebsiella spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Pseudomonas spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Enterobacter spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Serratia spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Hafnia spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	
Enterococcus spp	 < 1,0 x 10 ⁴	CFU/g faeces	10 ⁶ - 10 ⁷	



Determination of anaerobic bacteria

Name	Your value	Unit	Reference value	Scale
Bifidobacterium spp	 3,0 x 10 ⁹	CFU/g faeces	10 ⁹ - 10 ¹¹	
Bacteroides spp	 1,0 x 10 ⁹	CFU/g faeces	10 ⁹ - 10 ¹¹	
Lactobacillus spp	 < 1,0 x 10 ⁴	CFU/g faeces	10 ⁵ - 10 ⁷	
Clostridium spp	 < 1,0 x 10 ⁵	CFU/g faeces	< 1,0 x 10 ⁵	









Mycological stool examination

Name	Your value	Unit	Reference value	Scale
Candida spp	 < 1,0 x 10 ³	CFU/g faeces	< 1,0 x 10 ³	
Candida albicans	 < 1,0 x 10 ³	CFU/g faeces	< 1,0 x 10 ³	
Yeast	Negative	Negative		
Geotrichum candidum	 < 1,0 x 10 ³	CFU/g faeces	< 1,0 x 10 ³	





pH value

Name	Your value	Unit	Reference value	Scale
pH	 6.87		5,8 - 6,5	





Digestive residues

Name	Your value	Unit	Reference value	Scale
Quantitative determination of fat	 5.41	g/100g	< 3,5	
Quantitative determination of nitrogen	 0.80	g/100g	< 1,0	
Quantitative determination of sugar	 3.02	g/100g	< 2,5	
Quantitative determination of water	 74.36	g/100g	75 - 85	



Detection of indigestion


Name	Your value	Unit	Reference value	Scale
Pancreatic elastase	 396.71	µg/g	> 200	
Bile acids in stool	 79.54	µmol/l	< 70	

Detection of malabsorption

Name	Your value	Unit	Reference value	Scale
Calprotectin	 30.74	mg/l	< 50	
Alfa-1 antitrypsin	 31.54	mg/dl	< 27,5	

Mucosal immunity

Name	Your value	Unit	Reference value	Scale
Secretory IgA	 167.00	CFU/g faeces	510 - 2040	





Extended Information

Explanation of your test results

The intestinal flora includes the bacteria and other microorganisms that are naturally present in the intestine. Our intestinal flora has a broad influence on human health and disease. It modulates the immune system, provides the organism with vitamins, participates in the digestion of food components, supplies intestinal epithelium with energy by producing short-chain fatty acids and stimulates intestinal peristalsis.

Determination of aerobic bacteria

Escherichia coli



Name	Your value	Unit	Reference value	Scale
Escherichia coli	 8,0 x 10 ⁵	CFU/g faeces	10 ⁶ - 10 ⁷	

Escherichia coli belongs to the putrefactive bacteria. Putrefactive bacteria are bacteria that thrive in an alkaline environment and produce ammonia. At high pH, they grow and can compete with other good bacteria.

Putrefactive bacteria grow by metabolizing protein and fat, which produces toxic metabolites such as ammonia, indole, skatole and hydrogen sulfide. This can damage the intestinal mucosa and may lead to an increase in the pH of the colon (> 8.5) in the long term.

Elevated levels are often caused by increased activity in the immune system. If the immune system is characterized by an increased amount of putrefactive bacteria, attempts to reduce the pH value in the intestine (e.g. by supplementing with probiotics) can promote the reconstitution of the immune system and inhibit the growth of putrefactive bacteria. Less toxic metabolites are produced, which also reduces the damage on the intestine and other vital organs such as the liver and kidneys.

Escherichia coli Biovare

Name	Your value	Unit	Reference value	Scale
Escherichia coli Biovare	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	

Escherichia coli Biovare are subspecies of Escherichia coli that does not normally occur in the intestine. There are several types such as hemolytic, mucosal and lactose negative E. coli. It is not uncommon to find hemolytic or mucosal E. coli while looking for verotoxin E. coli and a low pH.

Elevated levels can also be caused by increased mucosal density. This is mainly due to too low production of secretory IgA.

Proteus spp

Name	Your value	Unit	Reference value	Scale
Proteus spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	



Elevated levels of Proteus species often occur in Crohn's disease and urinary tract infections.

Klebsiella spp

Name	Your value	Unit	Reference value	Scale
Klebsiella spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	



High levels indicate an increased inflammatory activity in the intestine. Klebsiella can release histamine and may cause abdominal pain, bloating, gas and diarrhea. Elevated levels of Klebsiella is common after prolonged use of antibiotics.

Pseudomonas spp

Name	Your value	Unit	Reference value	Scale
Pseudomonas spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	

Elevated levels of Pseudomonas species may occur in inflammation of the intestine. Elevated levels can lead to diarrhea and loose stool.

Enterobacter spp



Name	Your value	Unit	Reference value	Scale
Enterobacter spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	

Enterobacter species belong to the facultative bacteria. Facultative bacteria are bacteria that thrive in an alkaline environment and prefer anaerobiosis. At high pH they grow and outcompete with other good bacteria. Facultative bacteria primarily metabolize protein and fat, which produces metabolites such as ammonia, indole, skatole and hydrogen sulfide. This can damage the intestinal mucosa and may lead to an increase in the pH of the colon (-4.5) in the long term.

Enterobacter species are common bacteria in polluted water and are often found in the intestinal tract of humans and animals. High levels are often associated with an inflammatory process in the intestine/mucosa. Enterobacter are opportunists that can occur as pathogens in hospitals where they cause infections in people with weakened immune systems.

If the intestine/lumen is characterized by an increased amount of facultative bacteria, one can try to reduce the pH value in the intestinal lumen by supplementing with probiotics. This promotes the recolonization of the intestinal acid/bacterium flora and inhibits the growth of facultative bacteria. Low levels metabolites (ammonia) are produced, which alleviates the damage on the intestine and other vital organs such as the liver and kidneys.

Serratia spp

Name	Your value	Unit	Reference value	Scale
Serratia spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	



Serratia species occur in a variety of collections. Elevated levels also often occur in people with diarrhea.

Hafnia spp

Name	Your value	Unit	Reference value	Scale
Hafnia spp	 < 1,0 x 10 ⁴	CFU/g faeces	< 1,0 x 10 ⁴	

Elevated levels of Hafnia species can occur in people with diarrhea and can occur (often in people with mechanical intestine systems).

Enterococcus spp



Name	Your value	Unit	Reference value	Scale
Enterococcus spp	 < 1,0 x 10 ⁴	CFU/g faeces	10 ⁶ - 10 ⁷	

Enterococcus produces wide metabolic products and antibodies/substances that prevent foreign bacteria colonization in the small intestine (parasitosis in small intestine/colonization).

Low levels promote the occurrence of endogenous infections by disrupting intestinal barrier function.

Determination of anaerobic bacteria

Bifidobacterium spp



Name	Your value	Unit	Reference value	Scale
Bifidobacterium spp	 3,0 x 10 ⁹	CFU/g faeces	10 ⁹ - 10 ¹¹	

Bifidobacteria are needed, among other things, to protect against potentially pathogenic microbes. For the synthesis of short-chain fatty acids and to stimulate IgA⁺ (gut-associated lymphoid tissue), as well as to break down complex carbohydrates and indigestible fiber. They are part of the so-called flora and mainly use carbohydrates as their fuel substance. They produce short-chain fatty acids (acetate and butyrate), which lower the pH value in the intestine/lumen and therefore inhibit the growth and spread of pathogenic bacteria.

High bifidobacteria levels is determination of colonization resistance in the colon. Bifidobacteria develop a microbial barrier against infection by disrupting mucosal receptors. These counteract the colonization and spread of pathogenic bacteria, yeast or parasites.

Low levels can be due to a mechanical intestine/mucosa or too low fiber intake. Excessive levels are rare. Bifidobacteria is a very common bacterial strain used in probiotics. Bifidobacteria is also well supported by prebiotics.

Bacteroides spp



Name	Your value	Unit	Reference value	Scale
Bacteroides spp	 1,0 x 10 ⁹	CFU/g faeces	10 ⁹ - 10 ¹¹	

Bacteroides species belong to the genus Bacteroides. They tend to break down indigestible fiber and produce butyrate/butyrate acid.

As with *Bacteroides*, most bacteroides have a dense colonization resistance in the colon. *Bacteroides* develops a microbial barrier against infection by occupying mucosal receptors. These counteract the colonization and spread of pathogenic bacteria, yeast or parasites.



The low levels of Bacteroides species indicates disturbance in the microbial ecosystem and can lead to overgrowth of other microbes. The low levels can also reduce the production of short-chain fatty acids. A diet rich in fiber and probiotics is recommended at low levels.

Lactobacillus spp

Name	Your value	Unit	Reference value	Scale
Lactobacillus spp	 < 1,0 x 10 ⁴	CFU/g faeces	10 ⁵ - 10 ⁷	

Lactobacillus is a bacterial species that occur in the small intestine and forms a beneficial ecosystem and metabolite products that prevent foreign bacteria from colonizing in the small intestine. The low levels can occur in a variety of other things, neurodegeneration, food allergies or food intolerances. Excessive levels can occur with impaired digestive capacity.

Clostridium spp


Name	Your value	Unit	Reference value	Scale
Clostridium spp	 < 1,0 x 10 ⁵	CFU/g faeces	< 1,0 x 10 ⁵	

Clostridium species contain over one hundred different sub-species. Most species are not virulent and can have positive effects on good intestinal health. Elevated levels may occur in people with autism, as well as in IBS. Even low levels can occur in people with IBS.

Mycological stool examination


The mycological stool examination is used to detect any overgrowth of yeast and fungi. A possible fungal overgrowth is a result of unfavorable conditions in the intestine, which in turn is often due to a weakening of the bacteria.

Candida spp

Name	Your value	Unit	Reference value	Scale
Candida spp	 < 1,0 x 10 ³	CFU/g faeces	< 1,0 x 10 ³	

Candida species are typically part of the commensal flora, but can become pathogenic in people with weakened immune systems and may cause vaginal fungal infections.

Candida albicans

Name	Your value	Unit	Reference value	Scale
Candida albicans	 < 1,0 x 10 ³	CFU/g faeces	< 1,0 x 10 ³	

Candida albicans belongs to the group of facultative pathogenic yeast which, under certain circumstances, intensively multiplies and can cause mycotic mycoses. Candida albicans accounts for 80-90% of all Candida mycoses.

Organisms that come in contact with human mucous membranes are often colonized with yeast. To prevent recurrence of infections, regular use of oral hygiene products should be practiced regularly.



If the intestine is characterized by an increased amount of putrefactive bacteria, one can try to reduce their value in the intestine/lumen by supplementing with pre- or probiotics. This promotes the reconstitution of the intestinal ecosystem and inhibits the growth of putrefactive bacteria. Lactobacilli metabolite byproducts are produced, which also reduce the damage on the intestine and other vital organs such as the liver and kidneys.

Yeast

Name	Your value	Reference value
Yeast	Negative	Negative

Yeast is a type of fungus that normally occurs in the digestive system, but can cause disease if they overgrow. Negative = no overgrowth. Positive = overgrowth. Your result should be negative.



Geotrichum candidum

Name	Your value	Unit	Reference value	Scale
Geotrichum candidum	 < 1,0 x 10 ³	CFU/g faeces	< 1,0 x 10 ³	

Geotrichum candidum is a fungus that can be isolated from soil, waste, unwashed vegetables, fruits and dairy products. This fungus is also often detected in water and stool.

Geotrichum can cause disease in the case of a weakened immune system, fungus in a child's treatment or immunosuppressive treatment. This affects the oral cavity, nose and throat in addition to the effect on the intestine.

The properties of the stool - pH-value



Name	Your value	Unit	Reference value	Scale
pH	 6.87		5,8 - 6,5	

The pH value of the stool can indicate if there are any conditions of stress, decay or fermentation in the intestine. A too low pH value often occurs in combination with a dysfunctional intestine. For example, sugars are metabolized to lactic acids, which can cause bacteria to overgrow in the stool. A too high pH value can be due to excessive amounts of protein, which can stimulate certain intestine bacteria to produce ammonia and other metabolic products, thus raising the pH value of the stool.

An acid content that is too effective (abnormal pH) can be stabilized by the addition of probiotics, in combination with a diet rich in fiber and low incidence in fat and protein.



Digestive residues

Quantitative determination of fat

Name	Your value	Unit	Reference value	Scale
Quantitative determination of fat	 5.41	g/100g	< 3,5	



Elevated levels of fat in the stool may be due to diet (high fat diet) or indicate disorders of fat breakdown or absorption. It may be due to disorders in bile secretion, bile acid absorption or incomplete secretion of lipolytic enzymes from the pancreas.

Quantitative determination of nitrogen

Name	Your value	Unit	Reference value	Scale
Quantitative determination of nitrogen	 0.80	g/100g	< 1,0	



Elevated levels may indicate a disturbance in the breakdown or absorption of protein in the small intestine.

Quantitative determination of sugar

Name	Your value	Unit	Reference value	Scale
Quantitative determination of sugar	 3.02	g/100g	< 2,5	

Elevated levels are often due to carbohydrate intolerance. Common causes are lactose malabsorption or lactose intolerance. Intolerance or carbohydrate malabsorption may also be due to secondary pancreatic insufficiency with reduced secretion of carbohydrate-digesting enzymes.



Quantitative determination of water

Name	Your value	Unit	Reference value	Scale
Quantitative determination of water	 74.36	g/100g	75 - 85	

Low levels indicate a delayed passage through the intestine or constipation. Elevated levels indicate an accelerated passage or diarrhea loose stools.



Detection of indigestion

Pancreatic elastase

Name	Your value	Unit	Reference value	Scale
Pancreatic elastase	 396.71	µg/g	> 200	

Pancreatic elastase is a protease enzyme that is secreted by the pancreas and breaks down elastin, the most common protein in elastin fibers, and digests other proteins such as fibrin, hemoglobin and albumin. Low levels occur with weakened pancreatic function, gallstones, cystic fibrosis and can occur in vegetarians.

Bile acids in stool



Name	Your value	Unit	Reference value	Scale
Bile acids in stool	 79.54	µmol/l	< 70	

Bile acids are responsible for, among other things, the separation of cholesterol, absorption of fatty acids and fat-soluble vitamins in the small intestine and stimulate bowel movements. Normally, approximately 80% of the bile acids are recycled for reuse (enterohepatic recirculation). A possible cause for low levels in the enterohepatic recirculation leads to an increased secretion of bile acids and leads to secretion of fluid and electrolytes in the small intestine in combination with an increased intestinal permeability and intestinal permeability leading to diarrhea.

Elevated levels indicate impaired ability of breaking down fat in the food.

Detection of malabsorption

Calprotectin



Name	Your value	Unit	Reference value	Scale
Calprotectin	 30.74	mg/l	< 50	

Calprotectin is a marker for inflammation in the gastrointestinal tract. Elevated levels may be due to certain drugs, diseases, food allergies, intestinal infections, celiac, inflammatory bowel disease, ulcers, diverticulitis or colorectal cancer. In case of prolonged elevated levels you should consult a doctor for a proper evaluation.

The calprotectin levels in the stool closely correlates with the activity and extent of inflammatory mucosal changes. Very high calprotectin levels are found, for example, in active chronic inflammatory gastrointestinal diseases, mucosal ulcers or extensive ulcerative colitis lesions. People with chronic diffuse gastrointestinal diseases, food allergies or malabsorption instead show lower or moderately elevated levels (50-150 mg/l).

Regardless of the cause of the elevated levels, one should try to achieve healing of the intestinal mucosa. According to studies, this is achieved particularly well by the addition of phosphatidylcholine (lecithin), which also supports the formation of an effective mucosal barrier by stabilizing and strengthening the mucosa. Because intestinal microbes can partially convert phosphatidylcholine to TMA (trimethylamine), which is further metabolized in the liver to TMAO (trimethylamine N-oxide) and can promote vascular inflammation and atherosclerosis, it is advisable to reduce the dose of lecithin and increase the desired effect by adding Alpha-1 antitrypsin, butyrate and L-glutamine or proбиоти.

Alfa-1 antitrypsin



Name	Your value	Unit	Reference value	Scale
Alfa-1 antitrypsin	 31.54	mg/dl	< 27,5	

Alpha-1 antitrypsin regulates inflammatory response by blocking enzymes released by leukocytes and macrophages.

Elevated levels of alpha-1 antitrypsin indicate an inflammatory irritation of the intestinal mucosa. This can lead to a decrease in the absorption of nutrients from food. There is often a link between elevated alpha-1 antitrypsin levels and an increased permeability of the intestinal mucosa, which in turn can lead to an increased load on the body's systemic defense system.

Regardless of the cause of the elevated levels, one should try to achieve healing of the intestinal mucosa. According to studies, this is achieved particularly well by the addition of phosphatidylcholine (lecithin), which also supports the formation of an effective mucosal barrier by stabilizing and strengthening the mucosa. Because intestinal microbes can partially convert phosphatidylcholine to TMA (trimethylamine), which is further metabolized in the liver to TMAO (trimethylamine N-oxide) and can promote vascular inflammation and atherosclerosis, it is advisable to reduce the dose of lecithin and increase the desired effect by adding Alpha-1 antitrypsin, butyrate and L-glutamine or proбиоти.

Mucosal immunity

Name	Your value	Unit	Reference value	Scale
Secretory IgA	 167.00	CFU/g faeces	510 - 2040	

Secretory IgA neutralizes viruses and prevents pathogens from attaching to the surface receptors in the intestinal mucosa. Secretory IgA thus belongs to our "first defense" against antigens and pathogens in the gastrointestinal tract and respiratory tract and protects against exposure to food-borne antigens.

Low levels of sIgA lead to an increased amount of antigens and that ongoing infections can no longer be managed, which leads to a chronic susceptibility to infections. This means that low levels are associated with diseases such as allergies, increased susceptibility to infection, immunosuppressive conditions or fungal infection in the gut. Low levels of sIgA may also occur in case of protein deficiency.

High levels may indicate an increased activity of defense reactions in the intestinal mucosa, which may be due to inflammation or allergic processes.

At low levels, supplements of probiotics are recommended.

How can you use the results

5R is a treatment protocol developed by the Institute for Functional Medicine in the United States to restore gastrointestinal health and other problems that originate therefrom. If you have any imbalances, we recommend that you follow the program as described below. The program usually takes between 3-6 months to complete.

1. Remove

Find and exclude any microbes, foods or toxins that may interfere with normal intestinal function and exclude them. Examples of antimicrobial herbs are Allicin/garlic extract, berberine, caprylic acid, grape seed extract, olive leaf extract and oil of oregano. If you suspect food intolerances or a toxic load, we offer food intolerance tests and heavy metals tests.

2. Replace

Support the digestive system with suitable digestive support. For example, betaine HCl with pepsin, digestive enzymes or herbs such as cardamom, cinnamon, fennel, gentian, ginger and turmeric.

3. Reintroduce


Balance the intestinal flora with dietary fiber, prebiotic foods and probiotic supplements. Choose probiotics according to the intestinal bacteria you are low on according to the test results.

4. Repair

Support the intestinal mucosa. Examples of substances to support the intestinal mucosa are aloe vera, B-vitamins, essential fatty acids, L-glutamine, Marshmallow root and vitamin A.

5. Rebalance - Remedy the lifestyle

Balance your lifestyle; diet, exercise, sleep and stress.



Other information

Leaky gut

The intestine is of great importance for the health of the body. In addition to absorbing nutrients, it has a very important function for our immune system - the intestinal mucosa acts as a protective barrier between the body and our environment.

In order to perform its tasks correctly, the intestines need:

- A balanced bacterial population (our intestinal flora/microbiome)
- Adequate production of substances that support the intestinal mucosa
- Adequate levels of secretory IgA (sIgA)
- An undamaged intestinal epithelial cell layer

The intestinal flora is important in fighting endogenous infections. The commensal intestinal flora is a barrier that makes it difficult for unwanted bacteria from outside of the body to settle in the intestine. This is called colonization resistance.

The mucus produced by the intestinal mucosa is important for the transport of food in the intestines. It protects the mucosa, is important for its integrity and is therefore the basis for maintaining normal permeability.

In addition, sIgA is transferred to and distributed in the intestines with the mucus. sIgA has the important job of binding antigens of all types in the intestinal lumen, for example bacteria, viruses, large macromolecules, etc. In this way it reduces the stress on the mucous membrane while keeping the antigen inside the intestines, prevents it from coming into contact with the mucous membrane and leads it out of the body.

sIgA thus controls the antigens in the body and the frequency of infection. It is therefore very important for the immune system.

The controlled (nutrient) substance passage from the intestinal lumen into the circulatory system takes place at the intestinal epithelial layer. First, this is done by the cells themselves (transcellular) and secondly between the cells (paracellular). The transcellular transport takes place without further disturbance of the cells via simple diffusion, via specialized receptors or by enclosing the substances by invaginating the cell membrane (endocytosis). The paracellular passage through the single epithelial layer, on the other hand, takes place at so-called "tight junctions". These are parts of a protein network, which encloses the cells and seals its gaps more or less. In some places, at these "tight junctions", the network allows a controlled passage of liquid and dissolved or colloidal substances.

If one or even several of the above conditions for a healthy gastrointestinal system are missing, various health problems will soon arise. If the intestinal flora is not in good condition, for example, pathogens can spread more easily and cause infections. If sIgA is not available in sufficient amounts, the antigenic load on the immune system increases. If the epithelial layer and its tight junctions do not hold tight, large amounts of unwanted substances can pass into the body. These issues can spread over time and cause problems in different parts of the body – especially autoimmune diseases.

Leaky gut and its significance

A well-functioning absorption capacity in the intestine is crucial for providing the body with nutrition. On the other hand, the intestinal mucosa protects the body from pathogenic bacteria, contaminants, etc. Therefore, a controlled permeability in the intestines is of great importance for our health.

If the intestinal permeability increases, larger amounts of substances pass into the circulatory system. The increased amount of substances can have a negative effect on the body in the long-term. The consequence of this is a reaction of the immune system against these substances: first, the intestinal mucosa is infected, which in the long run leads to damage to the intestinal mucosa, which further increases the intestinal permeability and a vicious cycle is initiated.

In addition, there are other immunological reactions: if the immune system fights against harmless food components, which pass into the body due to the "leaky gut", this can lead to food allergies or food intolerances and more problems will occur over time. People with a "leaky gut" can - in the long run - produce antibodies, which link to the body's own organ surfaces and suddenly the body begins to fight against itself. Simply explained, this means that a permeable intestine can cause autoimmune diseases and this has been proven in studies for, among other things, type 1 diabetes, multiple sclerosis and rheumatoid arthritis. One can only assume that a permeable intestine can be a (contributing) cause of many more diseases.

Underlying causes of increased intestinal permeability

First of all, various intestinal diseases affect the permeability of the intestine. Diseases such as Crohn's and ulcerative colitis always lead to increased intestinal permeability. Celiac disease and various types of intolerances such as lactose, fructose, etc. also affect intestinal permeability. Exocrine pancreatic insufficiency (lack of enzymes to break down food) also damages the intestinal mucosa and increases its permeability. Furthermore, infections, imbalances in the intestinal flora, toxins, as well as acute and chronic physical and mental stress can affect the permeability of the intestinal mucosa due to the increase in stress hormones.

Causes of increased intestinal permeability can be:

- Inflammatory bowel disease
- Celiac disease
- Food intolerances
- Impaired pancreatic function
- Low secretory IgA
- Mental and physical stress
- Infections and incorrect intestinal colonization (bacteria, parasites, viruses and yeast)
- Alcohol
- Drugs
- Heavy metals
- Oxidative stress

Milder intestinal inflammation

Mild intestinal inflammation can be a consequence of different disease states and can occur in different parts of the intestine. Most often, the lower part of the small intestine or the upper part of the large intestine is affected, which limits nutrient uptake. Undigested food has a tendency to remain in the intestine, which can cause irritation.

In case of irritation in the stomach/intestine, it is good to take into account and exclude foods that you are intolerant to. A well-balanced diet should be developed to reduce the strain on the gut. This can be accomplished by eating several smaller meals a day instead of large, heavy meals, as well as avoiding certain foods that can sometimes interfere with proper gastrointestinal function:

- Gas-producing vegetables: Cabbage, leeks, onions, peppers, mushrooms, legumes
- Hot spices: Chili, pepper, paprika powder
- Foods high in fat: Fried food, fatty fish, lard, cheese, mayonnaise
- Foods high in sugar: Sweets, products containing natural and artificial sweeteners
- Beverage: Alcohol, coffee, carbonated beverages, and drinks that are too hot or cold

It can be beneficial to test for food intolerances and exclude foods that may further irritate. Different cooking methods can also facilitate or make the load on the intestine more difficult.

During asymptomatic periods, a diet rich in fiber is recommended to reduce inflammation. If you change your diet, remember to do it gradually.

Foods that are usually well tolerated:

- Amaranth, buckwheat, brown rice, millet and quinoa
- Eggs, fish, lean meat and poultry
- Cooked vegetables
- Potatoes and other tubers
- Porridge
- Ripe fruits

Nutritional deficiencies

Imbalance in the intestine can lead to impaired nutrient uptake. If you have any digestive issues, it can be a good idea to test yourself for any deficiency of vitamin B12, iron, magnesium and zinc.

MCT oil

MCT oil is composed of medium chain triglycerides, which are used in disorders of fat absorption and can be used as a readily available energy source. They can be more easily absorbed and digested independently of bile acids and digestive enzymes. Adding these fats should be done gradually to allow the intestines to adapt.

Fiber-rich food

Dietary fiber is a substance that the body can only use to a certain extent. Normally, fibers are excreted undigested, but can also be processed by the intestinal flora and provide beneficial metabolic products such as butyrate that are very valuable for the body.

- Fibers have a structure that requires them to be chewed more and longer. Chewing leads to the release of saliva, which in turn counteracts caries and tooth decay. Chewing a lot also leads to a faster feeling of satiety.
- Dietary fiber swells and binds water. This prolongs the feeling of satiety. It also increases the weight of the stool and stretches the intestinal walls, which triggers the muscles in the intestine to work and stimulates intestinal emptying. Improved bowel movements make it easier for toxic substances to be excreted from the body and reduce the risk of diseases linked to excessive toxic exposure.
- Fibers do not just bind water in the intestine. Bile acids produced from cholesterol are also bound up and excreted in the stool. The body is thus forced to form new bile acids from cholesterol and the cholesterol level is lowered.
- Dietary fiber is important for people with diabetes, as it has a positive effect on blood sugar.
- Dietary fiber stimulates the growth of beneficial bacteria and strengthens our natural intestinal flora while reducing the growth of unfavorable bacteria and harmful metabolic substances.

The general dietary recommendation for fiber intake is 25-35 grams of fiber per day for adults. In the table below you can see how much fiber different foods contain.

Type of food	Food	Fibers in g/100 g
Legumes (uncooked)	Broad beans	16,4
	Brown beans	16,4
	Mung beans	16,0
	Soybeans	15,3
	White beans	15,8
Fruits	Pomegranate	10,0
	Passion fruit	15,9
	Raisins	9,7
	Dried apricots	12,0
	Dried figs	18,5
Vegetables	Avocado	4,8
	Brussels sprouts	4,5
	Artichoke	5,0
	Sun-dried tomatoes	12,7
	Wheat sprouts	14,0
Grains	Oatmeal	15,0
	Oat bran	18,0
	Hard bread (wholegrain)	24,0
	Wheat bran	37,5
	Wheat germ	30,5

Support your intestinal flora through the right diet and beverage

We have a large number of bacteria in our gut. They support our immune system and intestinal mucosa by producing important metabolic products and other nutrients. Which substances the bacteria produce depends a lot on the substrate in the intestine, i.e. food components that can not be broken down by the body. Based on their metabolic products, they support either the acidification or the putrefaction process in the intestinal flora. If the relationship between these two groups is in balance, the intestinal environment is slightly acidic. Lactic acid cultures thrive best in acidic environments while putrefactive bacteria prefer an alkaline environment.

Undigested food residues after consuming high protein and fatty foods, are used by the putrefactive bacteria in the intestinal flora. The metabolic products produced by a dominant putrefactive intestinal flora can stress the liver, cause bloating and gas and in part have carcinogenic effects in the long term. Fiber-rich diets, on the other hand, support an acidic intestinal flora, which supports peristalsis and protects against pathogens as well as carcinogenesis.

A well-balanced diet can largely balance the relationship between acidification and putrefaction. The intake of fat and protein should therefore be kept at a balanced and moderate healthy level.

Some ways to reduce fat and protein intake

- Eat lots of fruit and vegetables
- Avoid excessive intake of high-fat foods
- Eat lean fish, red meat and poultry
- Replace white flour products with whole grains
- Prepare the food carefully by, for example, steaming or baking in the oven rather than broiling, grilling or frying
- Make sure the fats you eat are good fats

This test does not replace a medical consultation. Always search for healthcare if you experience severe symptoms.

